



JOBST# UltraSheer Vascular Support Hosiery

Features & Benefits

Easy to wear

- Seamless circular knitted for a soft, silky look and comfortable feel
- Reciprocated heel and toe reduces friction and offers better fit and wear
- Nylon and spandex yarns combined with a silicone blend softener for easier donning and comfort

Wide choice

- Maternity styles provide a little "extra" for the mother-to-be
- Available in Knee High, Thigh High, Pantyhose and Maternity styles
- Large range of colours including patterned

KNEE HIGH 15-20mmHg

75165-04	Black (Small)
75166-04	Black (Medium)
75167-04	Black (Large)
75168-04	Black (X-Large)
75165-25	Natural (Small)
75166-25	Natural (Medium)
75167-25	Natural (Large)
75168-25	Natural (X-Large)
75165-26	Suntan (Small)
75166-26	Suntan (Medium)
75167-26	Suntan (Large)
75168-26	Suntan (X-Large)
73633-00	Pattern (Small)
73633-01	Pattern (Medium)
73633-02	Pattern (Large)
73633-03	Pattern (X-Large)

THIGH HIGH 15-20mmHg

75169-04	Black (Small)
75170-04	Black (Medium)
75171-04	Black (Large)
75172-04	Black (X-Large)
75169-22	Natural (Small)
75170-22	Natural (Medium)
75171-22	Natural (Large)
75172-22	Natural (X-Large)
75169-23	Suntan (Small)
75170-23	Suntan (Medium)
75171-23	Suntan (Large)
75172-23	Suntan (X-Large)
73635-00	Pattern (Small)
73635-01	Pattern (Medium)
73635-02	Pattern (Large)
73635-03	Pattern (X-Large)

WAIST HIGH 15-20mmHg

75173-04	Black (Small)
75174-04	Black (Medium)
75175-04	Black (Large)
75176-04	Black (X-Large)
75173-12	Natural (Small)
75174-12	Natural (Medium)
75175-12	Natural (Large)
75176-12	Natural (X-Large)
75173-13	Suntan (Small)
75174-13	Suntan (Medium)
75175-13	Suntan (Large)
75176-13	Suntan (X-Large)
73637-00	Pattern (Small)
73637-01	Pattern (Medium)
73637-02	Pattern (Large)
73637-03	Pattern (X-Large)

MATERNITY 15-20mmHg

75177-01	Black (Small)
75178-01	Black (Medium)
75179-01	Black (Large)
75180-01	Black (X-Large)
75177-04	Natural (Small)
75178-04	Natural (Medium)
75179-04	Natural (Large)
75180-04	Natural (X-Large)

MATERNITY 20-30mmHg

75213-03	Black (Small)
75214-03	Black (Medium)
75215-03	Black (Large)
75213-06	Natural (Small)
75214-06	Natural (Medium)
75215-06	Natural (Large)

KNEE HIGH 20-30mmHg

75201-04	Black (Small)
75202-04	Black (Medium)
75203-04	Black (Large)
75204-04	Black (X-Large)
75201-19	Natural (Small)
75202-19	Natural (Medium)
75203-19	Natural (Large)
75204-19	Natural (X-Large)
75201-20	Suntan (Small)
75202-20	Suntan (Medium)
75203-20	Suntan (Large)
75204-20	Suntan (X-Large)
73634-00	Pattern (Small)
73634-01	Pattern (Medium)
73634-02	Pattern (Large)
73634-03	Pattern (X-Large)

THIGH HIGH 20-30mmHg

75205-04	Black (Small)
75206-04	Black (Medium)
75207-04	Black (Large)
75208-04	Black (X-Large)
75205-23	Natural (Small)
75206-23	Natural (Medium)
75207-23	Natural (Large)
75208-23	Natural (X-Large)
75205-24	Suntan (Small)
75206-24	Suntan (Medium)
75207-24	Suntan (Large)
73636-00	Pattern (Small)
73636-01	Pattern (Medium)
73636-02	Pattern (Large)
73636-03	Pattern (X-Large)

WAIST HIGH 20-30mmHg

75209-04	Black (Small)
75210-04	Black (Medium)
75211-04	Black (Large)
75212-04	Black (X-Large)
75209-13	Natural (Small)
75210-13	Natural (Medium)
75211-13	Natural (Large)
75212-13	Natural (X-Large)
75209-14	Suntan (Small)
75210-14	Suntan (Medium)
75211-14	Suntan (Large)
75212-14	Suntan (X-Large)
73638-00	Pattern (Small)
73638-01	Pattern (Medium)
73638-02	Pattern (Large)
73638-03	Pattern (X-Large)



KNEE HIGH 30-40mmHg

75241-04	Black (Small)
75242-04	Black (Medium)
75243-04	Black (Large)
75244-04	Black (X-Large)
75241-21	Natural (Small)
75242-21	Natural (Medium)
75243-21	Natural (Large)
75244-21	Natural (X-Large)

THIGH HIGH 30-40mmHg

75245-04	Black (Small)
75246-04	Black (Medium)
75247-04	Black (Large)
75248-04	Black (X-Large)
75245-20	Natural (Small)
75246-20	Natural (Medium)
75247-20	Natural (Large)
75248-20	Natural (X-Large)

WAIST HIGH 30-40mmHg

75249-04	Black (Small)
75250-04	Black (Medium)
75251-04	Black (Large)
75252-04	Black (X-Large)
75249-13	Natural (Small)
75250-13	Natural (Medium)
75251-13	Natural (Large)
75252-13	Natural (X-Large)

JOBST Relief

JOBST# Relief delivers the certainty of graduated sustained compression therapy.

Comfortable

- A generous toe allowance
- Soft silk texture
- Beige colour reduces visibility
- Available in extra-full calf sizing
- Open and closed toe designs to suit all user preferences
- Silicone blend softener makes Relief easy to don
- Wide comfortable support band that stays in place

Durable

- Reinforced heel
- Durable silicone construction

Safe/Effective

- Latex free construction reduces skin irritation
- Hosiery design produces graduated sustained compression



KNEE HIGH CLOSED TOE - Beige

114620A	20-30mmHg (Small)
114621A	20-30mmHg (Medium)
114622A	20-30mmHg (Large)
114698A	20-30mmHg (Large - Full Calf)
114623A	20-30mmHg (X-Large)
114624A	20-30mmHg (X-Large - Full Calf)
114630A	30-40mmHg (Small)
114631A	30-40mmHg (Medium)
114632A	30-40mmHg (Large)
114699A	30-40mmHg (Large - Full Calf)
114633A	30-40mmHg (X-Large)
114634A	30-40mmHg (X-Large - Full Calf)

THIGH HIGH OPEN TOE - Beige

78042-08	20-30mmHg (Small)
78043-08	20-30mmHg (Medium)
78044-08	20-30mmHg (Large)
78045-08	20-30mmHg (X-Large)
78059-08	30-40mmHg (Small)
78047-08	30-40mmHg (Medium)
78048-08	30-40mmHg (Large)
78049-08	30-40mmHg (X-Large)

THIGH HIGH CLOSED TOE - Beige

78042-09	20-30mmHg (Small)
78043-09	20-30mmHg (Medium)
78044-09	20-30mmHg (Large)
78045-09	20-30mmHg (X-Large)
78059-09	30-40mmHg (Small)
78047-09	30-40mmHg (Medium)
78048-09	30-40mmHg (Large)
78049-09	30-40mmHg (X-Large)

WAIST HIGH CLOSED TOE - Beige

114656A	20-30mmHg (Small)
114657A	20-30mmHg (Medium)
114658A	20-30mmHg (Large)
114659A	20-30mmHg (X-Large)

KNEE HIGH OPEN TOE - Beige

114625A	20-30mmHg (Small)
114626A	20-30mmHg (Medium)
114627A	20-30mmHg (Large)
114696A	20-30mmHg (Large - Full Calf)
114628A	20-30mmHg (X-Large)
114629A	20-30mmHg (X-Large - Full Calf)
114635A	30-40mmHg (Small)
114636A	30-40mmHg (Medium)
114637A	30-40mmHg (Large)
114697A	30-40mmHg (Large - Full Calf)
114638A	30-40mmHg (X-Large)
114639A	30-40mmHg (X-Large - Full Calf)

WAIST HIGH OPEN TOE - Beige

114660A	20-30mmHg (Small)
114661A	20-30mmHg (Medium)
114662A	20-30mmHg (Large)
114663A	20-30mmHg (X-Large)
114668A	30-40mmHg (Small)
114669A	30-40mmHg (Medium)
114670A	30-40mmHg (Large)
114671A	30-40mmHg (X-Large)

CHAPS STYLE OPEN TOE - Beige Left Leg

114672A	20-30mmHg (Small)
114673A	20-30mmHg (Medium)
114674A	20-30mmHg (Large)
114675A	20-30mmHg (X-Large)
114784A	30-40mmHg (Small)
114785A	30-40mmHg (Medium)
114786A	30-40mmHg (Large)
114787A	30-40mmHg (X-Large)

CHAPS STYLE OPEN TOE - Beige Right Leg

114676A	20-30mmHg (Small)
114677A	20-30mmHg (Medium)
114678A	20-30mmHg (Large)
114679A	20-30mmHg (X-Large)
114788A	30-40mmHg (Small)
114789A	30-40mmHg (Medium)
114790A	30-40mmHg (Large)
114791A	30-40mmHg (X-Large)

JOBST# Measuring

Sizing for medical compression stockings is based on the circumferences of the leg. It is best to measure early in the day before swelling occurs.

What is your sizing for Jobst medical compression stockings? Jobst UltraSheer, and Relief stockings, 15-20 (except for Men), 20-30, and 30-40 mmHg are sized by leg circumferences using the chart below. Sizing for men and women is the same.

All dimensions are in centimeters (cm).

Knee high stockings require an ankle and calf circumference. Thigh high stockings require ankle, calf and thigh circumferences. Waist, maternity and chaps styles require ankle, thigh and hip circumferences. The Large-Tall is for shoe size >12. Refer to each specific product for available

Size	**	Ankle	Calf	Thigh	Hip
Small	**	18-21	28-38	40-62	71-117
Medium	**	21-25	30-42	46-70	76-127
Large	**	25-29	32-46	54-78	81-137
X-Large	**	29-33	34-50	60-81	102-166
Large-Tall	**	25-29	32-46		
Large Full Calf	**	25-30	46-61		
X-Large Full Calf	**	30-36	46-61		