



Etac has made some fresh changes to the Clean product family – with a focus on making Clean products more modern and more hygienic.

The new Clean models will roll out from March 2017, once current 2016 stock has been exhausted.

Clean accessories are also being rationalised, with Green accessories now discontinued and only Grey accessories being made

Please contact your supplier to discuss stock levels and availability for specific products or accessories..

Fresh, new colours



etac® A cleaner Clean

Hygiene Enhancements

Cleaner design on castor attachments; no visible screws for better infection control and cleaning



Labels replaced by 'pad' printing
Reduced friction pattern on armrests for easier cleaning
Integrated footrest support in plastic castor bracket; less metal parts



Soft backrest accessory with buckle attachment to replace Velcro



Simplified Accessories

All Green accessories discontinued; only Grey accessories now available (to suit all colour models)



Rear wheel adaptor kit has white frame to fit both colours



Spare Parts – Castors include Attachments

Replacement Castors include attachments, and are backward compatible with old models

Current castors cannot be removed and remounted; when replacing castors, attachment must also be replaced.



BA3025	Etac Clean Mobile Shower Commode 55 cm, 4 Lockable Castors, Grey 15 kg (excludes Pan and Pan Holder) SWL 130kgs
BA3030	Etac Clean Mobile Shower Commode 55mm 4 Lockable Wheels Lagoon Green (excludes Pan and Pan Holder) SWL 130 kgs