



Access Health

MOTOmed Layson LA Supine Trainer

MOTOMED LAYSON LA SUPINE TRAINER

The MOTOMed Layson LA is a versatile medical device that provides patients in intensive care units, critical care and early rehabilitation with passive, motor-assisted and active training of the legs as well as the arms and upper body.

S.W.L
135 kg



MOTOmed®



EQUIPMENT

- The Layson LA can be easily and flexibly moved to the patient's bed and positioned in a space-saving manner due to its freely rotatable castors and the help of the optionally available expandable chassis
- The height adjustment and intuitive handling makes it is easily and effortlessly adaptable to the individual needs of each patient
- Large color touch display makes it easy to set programs and training parameters
- Plastic-coated safety foot shells with velcro
- QuickFix foot fastening system (Professional Model)
- TrainCare with plastic-coated leg guides and disinfectable fastening straps
- Knee bending adjustment incl. adjustment wheel
- Plastic-coated handles with quick release system
- Knee bending adjustment incl. adjustment wheel
- Removable operating panel (Professional Model)
- Patient remote control
- Hand fixation by wrist cuff (piece)

SPECIFICATIONS

- Total Length: 151 - 178cm
- Total Width: 62 -100cm
- Total Height: 118 - 166cm
- Unit Weight: 92kg

ITEM CODE

AP2162 MotoMed Layson LA Supine Trainer

AP2165 MotoMed Layson LA Supine Trainer - Professional

Ph: 03 9878 1399

| sales@accesshealth.com.au

| www.accesshealth.com.au



Quality
Approved
ISO 9001