



Access Health

Access General Ward Hospital Bed

GENERAL WARD HOSPITAL BED

The BAE502 series Electric Hospital Bed is a high quality adjustable electric hospital bed at an affordable price.




KEY FEATURES

- Battery backup
- Manual CPR handle and CPR buttons on the side rail control panel
- Caregiver patient side rails controls
- Caregiver foot end controls
- Lockout control on foot end
- 12° Trendelenburg and reverse Trendelenburg
- Head section and knee section can move at the same time for cardiac chair position
- 4-position IV pole mounts
- Blow-moulded headboard and footboard, single-step removal
- 4 centrally controlled castors (diameter: 125 mm)
- Standard length sleep deck – 1960 mm
- Drainage hook and IV pole
- Compatible with most:
 - standard mattresses
 - overbed tables in the low position
 - bedside cabinets



SPECIFICATIONS

- Bed length with standard sleep deck: 2300 mm
- Width (with resident assist rails up): 1050 mm
- Width (with resident assist rails down): 990 mm
- Height of sleep deck (in low position): 450 mm
- Height of sleep deck (in high position): 700 mm
- Sleep deck width: 900 mm
- Length: 1960 mm
- Max head elevation: 75°
- Max knee elevation: 40°
- Max Trendelenberg & Reverse Trendelenberg: 12°
- Castors: 125 mm



Easy control from side rails and bed end!

ITEM CODE

BE1505 ACCESS GENERAL WARD HOSPITAL BED

Ph: 03 9878 1399

sales@accesshealth.com.au

www.accesshealth.com.au

