




## SLEIGH BED

## ATTENDA<sup>®</sup> BEDROOM SERIES



**10 Year**  
Limited  
replacement  
warranty

## MATERIAL

- The slightly textured maintenance-free surface allows for easy cleaning. Rotationally moulded, specially formulated, fire-retardant, high-impact polyethylene with ultraviolet light stabilisers (to reduce fading). Chemically resistant to blood, vinegar, urine, faeces, salt solution and chlorine solution
- Premium TruColor™ material is fully compounded for superior colour, fire retardant and quality assurance. The unit is filled with rigid polyurethane foam for increased durability and sound absorption
- Made in the U.S.A

## OPTIONS

- Drainage holes
- Stainless steel plates for bolt-down application

## SPECIFICATIONS

- Size: 2100 x 800 x 720-950mm
- Backrest max upward angle: 90°
- Footrest max upward angle: 45°
- Trendelenburg: 18°
- Anti-trendelenburg: 18°

## DIMENSIONS

### BE2015

- Overall Length: 236cm
- Overall Width: 98cm
- Height: 64.5cm
- Surface Height: 33cm
- Weight: 56.5kg

### BE2020

- Overall Length: 236cm
- Overall Width: 98cm
- Height: 64.5cm
- Surface Height: 46cm
- Weight: 62.5kg

### ITEM CODE

BE2015 Attenda Sleigh Bed 33cm H (13") - Seclusion / Mental Health Bed

BE2020 Attenda Sleigh Bed 45cm H (18") - Seclusion / Mental Health Bed