



Access Health

LiftCare Floorline LC600 Electric Hospital Bed

LIFTCARE FLOORLINE LC600 ELECTRIC HOSPITAL BED



3 Year
Warranty
Other Parts

10 Year
Warranty
Frames

Easy to fold
& transport

KEY FEATURES

- Fully folding one-piece design
- 95mm-700mm 'true floor level' height range
- Electric Hi-Lo Backrest and Trendelenburg (& Reverse) standard
- Fully enclosed ABS plastic moulded covers
- 100mm castors
- 900mm wide 'smoothline' mattress platform
- Standard mattress platform 2000mm between head & footboards with optional 175mm Extension Kit & Bolster
- Ultra quiet operation using quality electrics
- Auto contour and lock out remote controls
- Colour: 'Vanilla' (Cream)

BENEFITS

- Best 'NO LIFT' practice; designed to prevent falls, minimise restraint, reduce risk of liability and litigation
- Improved infection control through Smoothline mattress platform
- ABS covers close 'seal' at top and smooth, clean lines on bed ends

SPECIFICATIONS

- Lowest mattress platform position: 95mm
- Highest mattress platform position: 700mm
- Sleeping length standard: 2000mm
- Sleeping length with Ext kit: 2175mm
- Mattress platform width: 900mm
- Overall bed weight: 100kgs
- Overall bed width: 905mm
- Overall bed length: 2150mm
- Extension kit overall bed length: 2325mm

ITEM CODE

BE3005 Liftcare Floorline LC600 Electric Hospital bed Floorline (See LTC)

Ph: 03 9878 1399 | sales@accesshealth.com.au | www.accesshealth.com.au

