

AIRCAST CRYO CUFF KNEE – ALL SIZES

The Knee Cryo/Cuff combines the therapeutic benefits of controlled compression to minimize hemarthrosis, oedema hematoma, swelling, and pain.



AIRCAST®

Completely
fit the knee
to provide
maximum
cryotherapy!

KEY FEATURES

- Anatomical cuff design for complete coverage of affected area
- Measured compression for patient comfort. Controlled cold eliminates the risk of tissue damage
- Ideal after for use post op. Measure at torso circumference (chest circumference) for sizing
- The cuff can be detached any time during treatment without interruption
- The Aircast Cryo/Cuff Cooler or IC Cooler Unit must be purchased separately



ITEM CODE

CC2030	Aircast Cryo Cuff - Medium Knee Fits 45cm - 60cm
CC2035	Aircast Cryo Cuff - Large Knee Fits 50cm - 80cm
CC2065	Aircast Cryo Cuff Knee Small Fits 25cm - 50cm