



Access Health

Kinetic Centura Shoulder CPM Machine

Kinetic Centura Shoulder CPM Machine is the ideal choice for hospital shoulder rehabilitation.

- Quality and number of movements
- Synchronisation system
- Flexible and evolutionary machine
- Soft chair for patient comfort

5
different
movements!

FEATURES

- Composed of 5 complementary and anatomical movements covering most of the shoulder pathologies
- Prevents joint stiffness, oedema, soft tissue contractures and muscles atrophy
- Patient set-up is anatomical, comfortable and quick with a colour-coded setup system
- Programmable hand control allows pre-programming of 16 different protocols
- Change for right and left is easy



SPECIFICATIONS

- Patient sizing: from 1.4 to 2m
- Maximum User Weight: 130kg
- Load reverse ensures patient protection when excessive force is exerted on the joint
- Adjustment of the inversion force: 6 levels
- Adjustment of the session time: from 1 minute to 24 hours
- Pause in movement limit: from 1 second to 15 minutes
- Adjustment of the speed from 50 to 140° per minute (5 levels)
- Dimensions: 56cm x 100cm x 76cm (variable height)
- Voltage: from 100 to 240 Volts 50/60Hz
- Electrical power: 50 VA
- Warranty: 2 years

ITEM CODE

CP2305 Kinetic Centura Shoulder CPM Machine