

OVERBED TABLE

PremiumLift antimicrobial table top slow the spread of microorganisms. Microorganisms include bacteria, viruses, protozoans, and fungi.

The thermoformed top is designed with infection control in mind. The surface is non-porous and gap free, which eliminates areas for water to collect and bacteria to grow.



FIXED



TILT



**3 Year
Warranty**

Ergonomic and Safe Auto-Lift Mechanism

PremiumLift Overbed Table unique design has an Auto-Lift mechanism allowing the patient to effortlessly raise the table top from any position. The table top will remain at the required position with no further adjustments and is lowered by using the handle with no bending required.

KEY FEATURES

- Ergonomic & Modern Design
The unique design allows the user to adjust the height with one hand no bending is required (Auto-Lift)
- Infection Control
The Antimicrobial table top slow the spread of microorganisms (Gap Free)
- Upgraded Castors
Heavy duty nylon castors (easy roll)
- Reduced Minimum Height
The minimum height has been reduced to 76cm
- ISO Certified & TGA Approved

SPECIFICATIONS

- Height Range: 760mm to 1030mm
- Castors: 50mm nylon - 2x lockable
- Weight: 15 kg
- Lift Mechanism: Ergonomic Auto Lift
- Split Tilt Top: 600 - 200 x 400mm
- Tilt Top Angle: 0 – 70 Degrees
- Non-Tilt Table Top: 800mm x 400mm



ITEM CODE

FR1126 Access Premium AutoLift Anti Microbial Overbed Table TILT Top Grey

FR1126F Access Premium AutoLift Anti Microbial Overbed Table FIXED Top Grey