

5 SECTION MANIPULATION COUCH

A popular selection for physiotherapists this five section couch offers strong smooth electric hi-lo operation through footswitch along with a generous range of motion for the head & backrest. This couch will allow correct patient positioning while reducing practitioner strain.

No Postural Drainage



Electric Postural Drainage operation by foot switch



KEY FEATURES

- Electric Hi Lo operation by foot switch
- Height adjustable padded armrests
- Head section has face hole and plug
- Adjustable feet for uneven surfaces
- 7.5cm retractable castors for safety and mobility
- Reliable TiMotion motor
- Supportive foam padding

CERTIFICATION

- CE, FDA, ISO13485-2016

UPHOLSTERY

- Heavy duty PVC vinyl that is fire, stain, mildew, oil, water and wear resistant

SPECIFICATIONS

- Total Length: 195cm
- Width: 66cm
- Headrest length: 43cm
- Mid section length: 51cm
- Backrest length: 96cm
- Adjustable backrest through gas operated strut to +70 Degrees
- Adjustable headrest through gas operated strut from -35 to +70 Degrees
- Height range from 50cm - 96cm in approx 30 sec
- Padding thickness: 5cm
- Foam Density: 36kg/m³
- Warranty: 5 Years Construction, 3 Years Electrics
- Weight Capacity: 250kg

Sleek design with rounded frame legs for safety!

ITEM CODE

- | | |
|--------|---|
| MC1010 | Access 5 Section Manipulation Couch - No Postural Drainage, No Lift Up Back, Up Down Arms |
| MC1011 | Access 5 Section Manipulation Couch - With Postural Drainage, Up Down Arm Rests |