

LINDA 3 PODIATRY CHAIR

Podiatric chair equipped with electronically controlled independent movements of height, backrest and leg-rest inclination. Perfect combination of elegant pure lines, reliability and ergonomics.

KEY FEATURES

- Electromechanical up/down movements, tilting backrest, adjustable leg rest
- Anatomical backrest
- Ergonomic PODO headrest
- Possible customization with different colors of polyurethane and fabric
- Fully left and right symmetry, for both left and right handed professionals
- Mobile high precision pedal control
- Emergency chair stops on backrest and leg rest
- Wrapping upholstery without seams
- Independent leg section movement
- RAL9002 standard chair painting
- Automatic return to "0" position after the treatment
- 3 programmable operating positions
- Minimum height 50cm, which allows easy access also for elder and disabled patients
- Wrapping upholstery without seams

SPECIFICATIONS

- Minimum height: 510mm
- Maximum height: 840mm
- Backrest tilting angle: -5/+70°
- Rotation angle: 325°
- Width: 660mm
- Length (without headrest): 15200mm
- Length (with headrest): 1900mm
- Gross weight: 145kg
- Net weight: 115kg
- Maximum lifting capability (patient + unit): 135+75
- Insulation classification: I
- Applied parts class: B

TECNODENT
ERGONOMIC INNOVATION

Ergonomic
and reliable!



OPTIONS

- Memory Foam or Hand Made Soft Foam upholstery
- Possibility to equip the chair with one/two rotating armrests
- Inox basin with support
- Different chair painting colors
- Magnetic leg rest pillow
- Children pillow
- Tray Holder with Tray
- Rotation



ITEM CODE

MC4110 Tecnodent Linda 3 Podiatry Chair Single leg support, foot control and Ergonomic Head Support