



Access Health

Monark 928G3 – Ergometer Upright Step Thru Bike

MONARK 928G3 – ERGOMETER UPRIGHT STEP THRU BIKE

The Monark 928G3 Ergometer Upright Step Thru Bike is designed for comfort and accessibility, featuring a step-through design for easy mounting and dismounting. It boasts a sturdy frame and adjustable seat for optimal posture and comfort during workouts. Equipped with a clear LCD display, it tracks essential metrics such as time, distance, and heart rate. Ideal for users seeking a reliable upright exercise bike, the 928G3 combines durability with user-friendly features to support effective fitness routines at home or in gym settings.

KEY FEATURES

- Super fast motor
- Easy entry frame
- Stand alone, 4 fixed test protocols
- Good cycle feeling
- Speed independent power control
- Easy to control and calibrate
- Remote control USB
- Easy to connect to PC Monark Test Software
- Workload range 15–700 W
- Comes with a Monark G3 Display-System
- The unique Monark G3 display unit can run 3 stand alone test protocols and 1 custom increment direct from the display interface
- Training is easy to control with the total summary of min and max data at the end
- Fully compatible with Monark Test Software.



S.W.L
250 kg



SPECIFICATIONS

- Resistance System: Spring
- Resistance Workload Range: 4–1, 400W
- Resistance Flywheel: 20 kg
- Applications: Heart Rate, Testing, Constant Power, Software, ECG, Mechanical Calibration
- Display: RPM, Heart Rate (ANT+ / 5 kHz Receiver), Time, Speed, Distance, Watt, METS
- Frame Type: Medical Walk-Through
- Handlebar: Monark Handlebar
- Seatpost and Saddle: Monark Standard
- Bottom Bracket: 68/122 mm
- Crank: Steel, 170 mm
- Pedals: 9/16", pedals with foot straps
- Max User Weight: 250 kg
- Length (mm): 1240 mm
- Width (mm): 500 mm at handlebar; 640 mm at support legs
- Height (mm): 1260 mm at handlebar; 780–1170 mm at saddle
- Unit Weight: 58 kg
- Power Requirements: AC 110V–220V | 50/60Hz

ITEM CODE

ITEM # MONARK 928G3 – ERGOMETER UPRIGHT STEP THRU BIKE

Ph: 03 9878 1399 | sales@accesshealth.com.au | www.accesshealth.com.au

