



Access Health

Dyaco 7.0R Rehabilitation Recumbent Bike

7.0R REHABILITATION RECUMBENT BIKE

The Rehabilitation Bike - 7.0R not only puts comfort at the forefront of a patient's experience, it also has been engineered for every single workout to be highly effective therapy for recovery.

Sitting in and exiting are made easy with its swivelling seat functionality. Perfect for deconditioned or assisted patients to be correctly aligned when sitting for their physical therapy session.

KEY FEATURES

- Adjustable Swivel Seat
- Adjustable Pedal Crank
- 5-750 Watt Range
- Low Step-Through Design
- 50 Resistance Levels
- 8 Position Swivel Seat
- 6 Position Recline Adjustments
- Download Workout Metrics

CONSOLE

- Display: LED & Dot Matrix window
- Programs: Manual, Hill Plateau, Interval profiles, VO2, HR, Custom, Isokinetic speed-based resistance (from 25 to 100 RPM), Symmetry
- Start Resistance: Yes (5 Watts)
- Hand Pulse: Yes
- Csafe: Yes

SPECIFICATIONS

- Seat adjustment: 8 position swivel seat adj., 6 position recline seat back, fore/aft positioning adj.
- Certificate: MDD 93/42/EEC Class Im cTUVus
- Max User Weight: 440lbs (200kg)
- Dimension: 57" x 30" x 51" / 145 x 77 x 130cm
- N.W: 180.1lbs / 81.7kg
- G.W: 198.9lbs / 90.4kg

Strengthen knees & ankles without the impact!

S.W.L
200kg

dyaco



Displays workout data



Customizable range of motion



ITEM CODE

RB2410 Dyaco Medical 7.0R Rehabilitation Recumbent Bike

Ph: 08 9873 1399

| sales@accesshealth.com.au

| www.accesshealth.com.au



Quality Approved
ISO 9001