



Access Health

Compact Wheelchair Scale

COMPACT WHEELCHAIR SCALE

The compact chair scale is specially designed for patients that have difficulty standing, allowing for the simple and precise weighing of patients.

S.W.L
250 kg



Optional
handrails can
be attached
for stability!

KEY FEATURES

- Collapsible folding frames for easy transportation
- Can be assembled easily by one person without tools
- The low-profile ramp leads to a non-slip platform designed especially for patients with difficulty standing
- By taring the weight of the wheelchair, a patient's exact weight can be accurately measured

SPECIFICATIONS

- 960 (H) x 1020 (D) x 590mm (W)
- Device weight: 17.6kg (38.8lb)

(NOTE: The device is NOT a wheelchair, and should not be pushed when a patient is sitting in the seat)



ITEM CODE

SS2046 Chair Scales Compact (MS5460) SWL 250kg

Ph: 03 9878 1399

| sales@accesshealth.com.au

| www.accesshealth.com.au



Quality
Approved
ISO 9001