




Access Health

TRIBUTE 7 Drawer Resuscitation Trolley - 6 Half & 1 Full Drawer

7 DRAWER RESUSCITATION TROLLEY - 6 HALF DRAWERS AND 1 FULL DRAWER

TRIBUTE

Double safety folds to all shelf edges!



Fill welding for infection control



Drawer dividers optional



125mm castors with 2 x locking castors

KEY FEATURES

- 6 x Half drawers 125mm deep, 1 lockable
- 1 x Drawer 250mm deep
- Full under bracing support for all shelves providing maximum stability and strength
- Double safety folds to all shelf edges
- Work surface with 3 rails
- 125mm castors, two locking
- 304 Grade stainless steel
- 125mm Satin finish

ITEM CODE

TM1250	TRIBUTE Stainless Steel 7 Drawer Resuscitation / IV set up Trolley 800 x 490 x H 1050mm - 6 Half Drawers and 1 Full Drawer, 1 x Lockable Drawer
TM1255	TRIBUTE Stainless Steel 7 Drawer Resuscitation / IV set up Trolley 900 x 490 x H 900mm - 6 Half Drawers and 1 Full Drawer, 1 x Lockable Drawer
TM1260	TRIBUTE Stainless Steel 7 Drawer Resuscitation / IV set up Trolley 1000 x 490 x H 900mm - 6 Half Drawers and 1 Full Drawer, 1 x Lockable Drawer

Ph: 03 9878 1399

| sales@accesshealth.com.au

| www.accesshealth.com.au

