

# Summit Anaesthesia Cart - 5 Drawer



## Access Health

### SUMMIT ANAESTHESIA CART - 5 DRAWER

A versatile anaesthesia trolley with central or seals drawer lock, soft closing drawers and smooth braked castors. Includes the extension table, storage rack, waste bin, steel basket, ABS drawer dividers & a non slip mat.



Durable and suitable for anaesthesia & many other applications!

### KEY FEATURES

- Extension Table
- Waste Bin
- 304 Stainless Steel Wire Basket
- ABS Divider
- Central Lock
- Castors with brakes
- Nonslip mat
- Luxurious push handle
- Soft closing drawers
- Drawer Sizes: 2 x 110mm, 1 x 420mm
- G.W. 35kg

### DIMENSIONS

- 650 x 480 x 980mm



#### ITEM CODE

TM8205 Summit Anaesthesia Cart - 5 Drawer - Includes Accessories

Ph: 03 9878 1399

| sales@accesshealth.com.au

| www.accesshealth.com.au



Quality  
Approved  
ISO 9001