



Access Health

Chattanooga Intellect Transport 2 Ultrasound Bundle

CHATTANOOGA INTELECT TRANSPORT 2 ULTRASOUND BUNDLE

The Intellect Transport 2 ULTRASOUND is an ultrasound therapy system used with or without an optional cart. Ultrasound frequencies of 1 or 3.3 MHz are available with pulsed or continuous duty cycles (10%, 20%, 50% and 100%).

The device comes standard with a 5cm² ultrasound applicator, with additional sizes are available. Coupling detection settings give the option to alert the clinician when the applicator has become decoupled from the patient by an audible beep and stopping the treatment timer.

KEY FEATURES

- **Large Color LCD Screen** – High resolution color screen gives real time information on the progress of the individual treatment. Clear, intuitive menus allow for quick set up and customisation of the device.
- **Portability** – Ergonomically designed handle, lightweight design and optional battery give you the versatility to treat patients wherever is required.
- **Ultrasound Applicator** – Improved coupling detection and monitoring for more effective and efficient treatments.
- Mobile, table top, or cart usage
- User-defined memory positions to save up to 15 custom protocols
- 1 and 3.3 MHz ultrasound with interchangeable applicator sizes
- Transmitters with LED coupling light in 2 cm², 5 cm², and 10 cm² sizes
- Quick connect magnetic cart attachment, allowing easy clinic to field use. No tools required for connecting or disconnecting



ULTRASOUND INDICATIONS

Application of therapeutic deep heat for the treatment of selected sub-chronic and chronic medical conditions such as:

Relief of pain, muscle spasms and joint contractures that may be associated with conditions such as:

- Adhesive capsulitis
- Bursitis with slight calcification
- Myositis
- Soft tissue injuries
- Shortened tendons due to past injuries and scar tissues
- Relief of sub-chronic and chronic pain and joint contractures resulting from:
 - Capsular tightness
 - Capsular scarring

ITEM CODE

US2206 CHATTANOOGA INTELECT TRANSPORT 2 ULTRASOUND BUNDLE

Ph: 03 9878 1399 | sales@accesshealth.com.au | www.accesshealth.com.au

