



Access Health

Access Wooden Walking Rails

WOODEN WALKING RAILS



FIXED



FOLDING

KEY FEATURES

- Walking Rails with shaped and lacquered wooden hand rails for user comfort
- Height adjustable from 70 - 100cm in 25mm increments
- Available in 3 or 4m lengths and in either folding/free standing or fixed styles
- Approx 75cm between rails on folding style

**CUSTOMISED
OPTIONS**



ITEM CODE

WR1350	Access Fixed Walking Rails 2m Wooden Hand Rail
WR1353	Access Fixed Walking Rails 2.5m Wooden Hand Rail
WR1356	Access Fixed Walking Rails 3m Wooden Hand Rail
WR1359	Access Fixed Walking Rails 4m Wooden Handrail

ITEM CODE

WR1403	Access Folding Walking Rails 2m Wooden Hand Rail
WR1406	Access Folding Walking Rails 2.5m Wooden Hand Rail
WR1409	Access Folding Walking Rails 3m Wooden Hand Rail
WR1412	Access Folding Walking Rails 4m Wooden Hand Rail

Ph: 03 9878 1399

| sales@accesshealth.com.au

| www.accesshealth.com.au



Quality
Approved
ISO 9001