



Access Health

Tecnonodent Linda Next ENT Chair

TECNODENT LINDA NEXT ENT CHAIR

Elegant and versatile, ergonomic and functional, Linda NEXT ENT is a premium chair intended for audiological, ENT and related diagnostic examinations.

STANDARD FEATURES

- Complete symmetry for right – or left-handed professionals
- Electromechanical lifting, lowering and backrest tilting
- Backrest movement synchronised with leg section
- Automatic extensible PU foot section
- Independent Trendelenburg movement
- Automatic built in programs for patient positioning: table, Trendelenburg, rinse, last and zero position
- 3 operator profiles with 3 programmable positions each
- Safety security STOPS on backrest, leg rest and controls
- SLOW MODE feature for precise movements
- 24 DC linear actuators
- Soft start and soft stop
- Joystick control
- Double-articulated headrest VALE
- Standard molded upholstery
- RAL 9003 painting
- Different PU and fabric colors

OPTIONS

- Independent leg section
- Right and/or left swiveling armrest
- Different headrests available according to specific treatment needs
- Seamless memory foam upholstery
- Handmade soft-foam upholstery with contrasting color stitching
- Backrest motorised height adjustment
- Wheels
- Rotation system around the vertical axis under the chair base (320°)
- Backup Battery
- IR remote control
- Children pillow
- Leg section foot protection
- Stainless steel lateral tool holder bars
- Articulated arm with large instrument tray
- Surgical armrests
- Small side instrument tray with armchair attachment
- Multi-articulated headrest model "LIGHT"
- COMFORT wide multi-articulated headrest
- Paper roll holder



ITEM CODE

MC4080 TECNODENT LINDA NEXT ENT CHAIR

Ph: 03 9878 1399 | sales@accesshealth.com.au | www.accesshealth.com.au

