



# 102

## Fixed Height Practice Couch

With a sturdy powder coated frame, two drawers and a paper roll all in the price it's hardly a surprise this practical couch is such a favourite. Our Otis range of trolleys are designed to complement the couch (with Otis 1 being designed to fit under the couch) to offer maximum storage. The addition of steps and a matching operators chair would form the ideal package for an examination room.

Features: 1 year unconditional guarantee (excludes upholstery), strong powder coated frame, paper roll holder, 2 smooth running drawers, failsafe backrest mechanism, maintenance free, fast delivery.

### STANDARD SPECIFICATION

- Robust powder coated frame.
- 2 smooth running drawers
- Paper roll holder.
- Failsafe manual backrest mechanism.
- Backrest adjustment from 0° to +85°.
- Adjustable foot for uneven floors.
- Fully washable/ fire retardant Ambla vinyl available in a range of colours.
- Couch goes horizontal for patient recovery in emergency situations.
- CE marked in conformity with MDD 93/42/EEC Annex IX Class1.

### STANDARD OPTIONS

- Otis 1 three drawer trolley (fits under couch).
- Steps.
- Operators chair.

### 102 DIMENSIONS

**Length:** 188 cm (74")  
**Width:** 63 cm (25")  
**Height:** 82 cm (32") approx



#### 102 Fixed height couch

Please specify motorised back rest or other optional extras if required.

### PACKED WEIGHTS AND DIMENSIONS FOR SHIPPING

Packed Weight: 30 Kgs  
 Shipping Dimensions: 190 x 66 x 35 cm

### DELIVERY PROMISE

Our normal delivery service is 7 – 10 days from receipt of order; however if you have a specific requirement we will always endeavour to meet it.

### FOR YOUR PEACE OF MIND

Our unique guarantee is for a LIFETIME on the mechanism and UNCONDITIONAL for 3 years on all components.. This even applies to accidental damage. We will not quibble. The only exclusion is damage to upholstery.

All of our products are CE marked. We are an ISO9001 audited company. Should you require any help concerning the most suitable product for your needs, we would be delighted to discuss it with you.

#### Plinth 2000 Ltd.

Wetheringsett Manor  
 Wetheringsett  
 Stowmarket  
 Suffolk IP14 5PP  
 United Kingdom

t: +44 (0)1449 767887  
 f: +44 (0)1449 766122  
 e: sales@plinth2000.com

[www.plinth2000.com](http://www.plinth2000.com)

**Plinth** | 2000



# VINYL

Upholstery

We use only the best quality vinyl proven over many years and **ALL** feature an anti microbial additive to help prevent cross infection.

All coverings are vinyl and washable, being fire retardant to BS 7176, BS 5852 part 2 (crib 5) and display excellent levels of wear. ALL include patented anti microbial technology proved effective in inhibiting the growth of E-Coli and Staphylococcus Aureus (MRSA), and all are suitable for cleaning and disinfecting with NaDCC based products.

Note: Unless specified BEIGE will be supplied.



Beige



Slate



Night

**Plinth 2000 Ltd.**

Wetheringsett Manor  
Wetheringsett  
Stowmarket  
Suffolk IP14 5PP  
United Kingdom

t: +44 (0)1449 767887  
f: +44 (0)1449 766122  
e: sales@plinth2000.com

[www.plinth2000.com](http://www.plinth2000.com)

**Plinth** | 2000