



40

Wide Neurology Couch

Developed in conjunction with the Bobath Centre in London, our 'double' plinth provides probably the most firm and stable treatment surface available. Being 102cm wide and lowering to just 41cm it is perfect for use in hospital gyms & physiotherapy departments. Now safely lifts 320kg from its lowest point.

Features: unique lifetime guarantee, easy electric operation, fully hoist accessible, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 320kg from its lowest height, full choice of accessories and colours.

STANDARD SPECIFICATION


- Either electric or hydraulic height adjustment available.
- Able to lift up to 320kg (50 stone) from its lowest height.
- Easily able to lift and support two people.
- Backrest adjustment by gas strut mechanism.
- Breathing hole and plug fitted as standard.
- Backrest angle adjustable from 0° to + 85° vertical.
- Hoist clearance of 14cm (5 1/2").
- Adjustable foot for uneven floors.
- Fully washable/ fire retardant Ambla vinyl available in a range of colours.
- Couch goes horizontal for patient recovery in emergency situations.
- Heavy duty retractable rubber tyre castors for excellent mobility.
- CE marked in conformity with MDD 93/42/EEC Annex IX Class 1.
- Electrical standard IEC601 part 1.

STANDARD OPTIONS

- Motorised backrest (recommended).
- Battery backup (electric version only).
- Cot sides.
- Paper roll holder.
- Other options available.

40 DIMENSIONS

Length: 189 cm (74")
 Width: 102 cm (40")
 Height: 41 cm (16") approx
 to 94 cm (37") approx



40 Electric
 40 Hydraulic

PACKED WEIGHTS AND DIMENSIONS FOR SHIPPING

Packed Weight: 130 Kgs
 Shipping Dimensions: 190 x 102 x 57 cm

DELIVERY PROMISE

Our normal delivery service is 7 – 10 days from receipt of order; however if you have a specific requirement we will always endeavour to meet it.

FOR YOUR PEACE OF MIND

Our unique guarantee is for a LIFETIME on the mechanism and UNCONDITIONAL for 3 years on all components.. This even applies to accidental damage. We will not quibble. The only exclusion is damage to upholstery.

All of our products are CE marked. We are an ISO9001 audited company. Should you require any help concerning the most suitable product for your needs, we would be delighted to discuss it with you.



Plinth 2000 Ltd.

Wetheringsett Manor
 Wetheringsett
 Stowmarket
 Suffolk IP14 5PP
 United Kingdom

t: +44 (0)1449 767887
 f: +44 (0)1449 766122
 e: sales@plinth2000.com

www.plinth2000.com

Plinth | 2000



VINYL

Upholstery

All coverings are vinyl and washable, being fire retardant to BS 7176, BS 5852 part 2 (crib 5) and display excellent levels of wear. Most include patented anti microbial technology proved effective in inhibiting the growth of E-Coli and Staphylococcus Aureus (MRSA), and all are suitable for cleaning and disinfecting with NaDCC based products.

Note: Unless specified WADI will be supplied.


STANDARD HARD WEARING LEATHER GRAIN - STAIN RESISTANT AND FAST DRYING



Snow



Brume



Nubian




Wadi



Fern




Spruce



Memphis




Cornflower



Rainforest

NEW COLOURS - BRIGHTEN UP OR STAND OUT



Apple Mint



Lilac




Candy



Cool Blue




Mulled Wine




Cocoa


NON SLIP - GRIPPY, MORE COMFORTABLE



Skye



Daviot



Brechin

Plinth 2000 Ltd.

Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk IP14 5PP
United Kingdom

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com

www.plinth2000.com

Plinth | 2000