



503T

Tilting Minor Surgery Couch

Based on the popular model 503, this tilting version offers all the same features with the additional benefit of a seat tilt which offers greater comfort for patients and better positioning for the practitioner. Ideally suited for Minor Procedures, General Examination and Chiropody use.

Features: motorised seat tilt, smooth electric operation, unique lifetime guarantee, easy operation, low running costs, high density fire retardant foam, hardwearing anti-microbial vinyl, able to lift 160kg from its lowest height, full choice of accessories and colours.

STANDARD SPECIFICATION


- Electric height and trendelenburg tilt adjustment as standard.
- Able to lift up to 160kg (25 stone) from its lowest height.
- Backrest adjustment by gas strut mechanism.
- Operated as standard by a hand switch (footswitch optional).
- Backrest & leg rest adjustment from -20° to +85° & seat tilt from 0° to +28°.
- Hoist clearance of 13cm (5").
- Adjustable foot for uneven floors.
- Fully washable/ fire retardant Ambla vinyl available in a range of colours.
- Couch goes horizontal for patient recovery in emergency situations.
- Robust retractable castors for excellent mobility and safety.
- CE marked in conformity with MDD 93/42/EEC Annex IX Class1.
- Electrical standard to IEC601 part 1.

STANDARD OPTIONS

- Foot switch.
- Arm tables.
- Paper roll holder.
- Breathing hole & plug.
- Cot sides.
- Motorised backrest.

503T DIMENSIONS

Length: 190 cm (73")
Width: 63 cm (25")
Height: 46 cm (18") approx
to 95 cm (37") approx



503T Electric

Please specify motorised back rest or other optional extras if required.

PACKED WEIGHTS AND DIMENSIONS FOR SHIPPING

Packed Weight: 95 Kgs
Shipping Dimensions: 190 x 66 x 64 cm

DELIVERY PROMISE

Our normal delivery service is 7 – 10 days from receipt of order; however if you have a specific requirement we will always endeavour to meet it.

FOR YOUR PEACE OF MIND

Our unique guarantee is for a LIFETIME on the mechanism and UNCONDITIONAL for 3 years on all components.. This even applies to accidental damage. We will not quibble. The only exclusion is damage to upholstery.

All of our products are CE marked. We are an ISO9001 audited company. Should you require any help concerning the most suitable product for your needs, we would be delighted to discuss it with you.

Plinth 2000 Ltd.

Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk IP14 5PP
United Kingdom

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com

www.plinth2000.com

Plinth | 2000



VINYL

Upholstery

We use only the best quality vinyl proven over many years and **ALL** feature an anti microbial additive to help prevent cross infection.

All coverings are vinyl and washable, being fire retardant to BS 7176, BS 5852 part 2 (crib 5) and display excellent levels of wear. ALL include patented anti microbial technology proved effective in inhibiting the growth of E-Coli and Staphylococcus Aureus (MRSA), and all are suitable for cleaning and disinfecting with NaDCC based products.

Note: Unless specified BEIGE will be supplied.



Beige




Grey




Black



Fawn



Chocolate




Brow




Damson




Wine



Ruby



Terracotta



Night




Ocean




Sky



Slate




Green



Sage



Lime



White

Plinth 2000 Ltd.

Wetheringsett Manor
Wetheringsett
Stowmarket
Suffolk IP14 5PP
United Kingdom

t: +44 (0)1449 767887
f: +44 (0)1449 766122
e: sales@plinth2000.com

www.plinth2000.com

Plinth | 2000