

# Plinth 2000

## Model 93B Phlebotomy Chair

### Overview:

Flexible phlebotomy chair with unrivalled patient comfort and access. The easy removal of the fully adjustable arm rests mean this chair can also be used for general examination. Optional accessories can also be added to expand the use of this chair.




### Standard Features

- UK Manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Hydraulic or electric operation
- Comfort foam
- Easy operation
- Low running costs
- Hardwearing anti-microbial vinyl
- Safe working load 250kg
- Retractable castors
- Phlebotomy arm rests
- Hoist clearance 130mm


### Overall Dimensions

Length: 177cm (70") approx  
Width: 63cm (25") approx  
Height: 52cm - 101cm approx  
(20" - 40")approx


### Upholstery Dimensions



### Recommended Options:




Stainless Tray




Battery Back-up




D Cushion



Motor Back



CPR Motor Back



Paper Roll Holder


NHS Approved Supplier



t: +44 (0)1449 767887  
e: sales@plinth2000.com  
www.plinth2000.com



## Operational Parameters



Model 93BE - Electric version

Model 93BH - Hydraulic version

## For Your Peace of Mind

- Lifetime frame & 5 year parts guarantee (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average lifespan of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2008 by LRQA

## Standards

- IP Rating - IPx4
- EEC Classification - MDD93/42/EEC Annex IX Class 1
- Electrical Standard - IEC601 Part 1 + 60601. P (EN)
- Vinyl Flame Retardant To CRIB 5 (BS7176 Including Foam)
- Foam To BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

## Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required

## Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note: Unless specified, beige will be supplied.



Plinth 2000 Ltd.  
 Wetheringsett Manor  
 Wetheringsett  
 Stowmarket  
 Suffolk, IP14 5PP

t: +44 (0)1449 767887  
 f: +44 (0)1449 766122  
 e: sales@plinth2000.com  
 www.plinth2000.com

**Plinth** | 2000

