

Plinth Patient Chairs



Model 503CDT (NV0923)

Strong and extremely well proven podiatry chair that has set the market standard. Offering excellent patient comfort, strength and stability, it also benefits from unrivalled accessibility for both patient and operator.

This model has a motor which tilts the whole chair backwards. This greatly enhances comfort and enables the foot end to rise to nearly 153 cm (60").

Features:

- Electric height adjustment
- Extending and divided legs
- All functions operated as standard by a hand switch (or optional footswitch)
- Back rest adjustable from -20° to +90°
- Stay level swing back arm rests for easy patient access - arms are also fully removable
- Hoist clearance of 130mm
- Adjustable foot for uneven floors
- Heavy duty retractable castors for excellent mobility
- Couch goes horizontal for patient recovery in emergency situations

Technical Data

Length:	1890mm
Width:	650mm
Height:	Adjustable: 460 to 950mm
Leg length:	Adjustable: to 200mm
Foot section:	Adjustable: 30° to +20°
Maximum load on seat:	160kg

Model 93CD (NV2400)

Designed to allow a virtual 90° drop on the leg rests whilst still maintaining complete stability and patient safety. This model has all the features of the 500 range but is stronger, more versatile and now features a smooth gas spring operated back and foot rests

Features:

- Strong and extremely well proven podiatry chair that has set the market standard.
- Offering excellent patient comfort, strength and stability, it also benefits from unrivalled accessibility for both patient and operator
- Electric height adjustment
- Stay level swing back arm rests for easy patient access (arms also fully removable)
- Hoist clearance of 130mm
- Adjustable foot for uneven floors
- Fully washable and fire retardant Ambla vinyl available in a range of colours
- Heavy duty retractable castors for excellent mobility
- Couch goes horizontal for patient recovery in emergency situations

Technical Data

Length:	1720mm
Width:	630mm
Height:	Adjustable: 460 to 960mm
Maximum load on seat:	225kg

