



Access Health

Physical Therapy Supplies

Unit 1, Rear
194-196 Whitehorse Road
Blackburn, VIC, 3130
ph: (03) 9878 1377
fax: (03) 9878 7099
sales@accesshealth.com.au

FORMA SPLINT PAN 1501



Model 1501

Two convenient sizes for a variety of applications

Both Model 1501 and Model 1502 splint pans feature a digital thermometer on the control panel. Hinged, insulated cover. Stainless steel tank is trimmed in tough black polycarbonate. All-brass drain system. Grounded hospital-grade plug and a high heat safety switch. Forma-Splint™ units are UL listed, C-UL listed, CE conforming and are warranted for one year against defective parts or manufacturing.

Model 1501 features a dual thermostat that is factory set at 175°F (high) and 160°F (low). Settings may be readjusted. Operating temperature is reached in only 90 minutes with 3" (7.6cm) of cold water—even faster with less. This full size thermal bath with spacious stainless steel splint pan accommodates splinting sheets up to 24" x 18" (61 x 46cm). Tank interior is constructed of stainless steel, which heats evenly and resists corrosion. The hinged and insulated stainless steel lid is functionally designed to open completely or one side only, with just one hand. The split lid helps decrease condensation and maintain constant water temperature. This design is ideal for hand splints, allowing access to half the pan.

Dual thermostatic controls can be adjusted individually to conform to particular requirements. Dual independent heating elements for both settings provide quick warm up and automatically convert to a more efficient circuit once the working temperature is reached.

Model 1501 specifications: Overall measurements are 72 x 56 x 18cm.
Inside dimensions are 62 x 47 x 13cm. Weight: 21kg.