



FORM FIT® KNEE HINGED LATERAL J

The Form Fit Knee Hinged Lateral J is designed to provide realignment of the patella. It offers support, stability and protection for users benefitting from realignment of the patella and/or added lateral support of the patella.

FEATURES AND BENEFITS

- Breathable material
- Removable hinges
- Easy pull-on
- Lateral J buttress
- Open popliteal

INDICATIONS

Patients that would benefit from realignment of the patella and/or patients that would benefit from an added lateral support of the patella, such as:

- Lateral patellar subluxation and dislocations
- Lateral patellofemoral malalignment
- Patellofemoral maltracking and pain
- Chondromalacia due to lateral maltracking

PART#	SIZE	MEASUREMENTS
7091X2	XSmall	Thigh Circ: 13-15.5" (33-39.4cm) Knee Circ: 12-13" (30.5-33cm) Calf Circ: 10-12" (25.4-30.5cm)
7091X4	Small	Thigh Circ: 15.5-18.5" (39.4-47cm) Knee Circ: 13-14" (33-35.6cm) Calf Circ: 12-14" (30.5-35.6cm)
7091X5	Medium	Thigh Circ: 18.5-21" (47-53.3cm) Knee Circ: 14-15" (35.6-38.1cm) Calf Circ: 14-16" (35.6-40.6cm)
7091X7	Large	Thigh Circ: 21-23.5" (53.3-59.7cm) Knee Circ: 15-17" (38.1-43.2cm) Calf Circ: 16-18" (40.6-45.7cm)
7091X8	XLarge	Thigh Circ: 23.5-26.5" (59.7-67.3cm) Knee Circ: 17-19" (43.2-48.3cm) Calf Circ: 18-20" (45.7-50.8cm)
7091X9	2XLarge	Thigh Circ: 26.5-29.5" (67.3-74.9cm) Knee Circ: 19-21" (48.3-53.3cm) Calf Circ: 20-22" (50.8-55.9cm)

X = 5 (Right Leg) or 6 (Left Leg)

Thigh Circumference: 6" above mid-patella

Calf Circumference: 6" below mid-patella