

Easy To Use Fall Prevention Products For Caregivers

Caregiver Alert Motion Sensing Exit Alarm



RP-TL-2700

Features:

- Chime Alert
- Adjustable alert volume
- Wide field of detection
- Use for bedsides, hallways and doorways

Portable Caregiver Alarm With Motion Sensor



RP-TL-2800

Features:

- Setup is quick and easy - insert batteries and it's ready to use
- Motion Sensor can be placed by bed or doorway
- Optional wall mount accessories included

Caregiver Pager With Bed Exit Alarm



RP-TL-5102MP

Features:

- Quiet - No alarming monitors
- Easy to set up - insert batteries and it's ready to use
- Pocket sized pager - Easy to carry
- Exit Alarm can be placed by bed or doorway to monitor residents activities

Caregiver Pager With Call Buttons



RP-TL-5102TP

Features:

- Audible alert pager with adjustable tone
- Lanyard strap included
- Quiet at the bedside

Low Cost Quiet, Wireless Fall Prevention Alarm



No Cords! Place monitor where you need it!

Low Cost Solution

The Economy CordLess® Exit Alarm provides affordable and effective cord-free wireless monitoring with a range up to 300 ft. Ideal for home or small facility use.

Easy To Use

The Economy CordLess® Exit Alarm is simple-to-use and works with all 433 CordLess® Sensor Pads & Mats, Nurse Call Buttons, and Motion Sensors. Monitor is light-weight and can be paired with up to six (6) wireless components to serve as a portable caregiver alert.

Remove Noise From The Room

Monitor can be mounted away from the bed, removing in-room alarm noise and reducing entanglements and tripping hazards. Adjustable volume and gentle chime alert won't startle residents and offers QUIETER fall prevention.

Key Features:

- Affordable & effective wireless monitoring
- Easy set up and operation
- Works with all 433 CordLess® Sensor Pads & Mats, Nurse Call Buttons, and Motion Sensors
- Monitor up to six (6) components with one alarm
- Adjustable volume – low, medium and high volume adjustment and gentle chime setting
- Battery Operated (3-C, *Not Included*)
- Optional AC Power Adapter (AC-04)



Economy CordLess® Exit Alarm
433-EC



Call or Email Us Today! 03 9878 1399 sales@accesshealth.com.au

Lowest Cost Bed & Chair Exit Alarm



Lowest cost fall prevention monitor on the market!

Full-Function / Great Value!

The Smart Caregiver TL-2100B Monitor is a Full-Function, Low Cost fall alarm monitor designed for use with all Smart Caregiver Corded Pads and Mats. This is a high quality product at a great, everyday value.

Auto-Reset Function

Monitor can be placed at the bedside or wheelchair, connected to a corded sensor pad. The sensor pad is placed under the resident. When the resident gets up and pressure is removed from the pad, the monitor will sound alerting the caregiver. When pressure is re-applied to the pad, the monitor will automatically silence and reset for use.

Monitor can be set to silent

Monitor can be set to silent when using an existing nurse call system and additional nurse call cord.



TL-2100B

Key Features:

- Adjustable volume control – Monitor can be set to silent if used with existing nurse call system.
- Auto-Reset: Monitor automatically resets when pressure is re-applied to connected pad or mat.
- Nurse call port – connects monitor directly to existing nurse call system.
- Low battery & status indicator lights.
- Battery Operated (1-9Volt, Not Included).
- Optional AC Power Adapter (AC-02).

