



# Access Health

Physical Therapy Supplies

Unit 1, Rear  
194-196 Whitehorse Road  
Blackburn, VIC, 3130  
ph: (03) 9878 1377  
fax: (03) 9878 7099  
sales@accesshealth.com.au

## Hayman Chair



Weight Capacity: 300kg

Height: Adjustable 430mm - 580mm in 25mm increments

Seat Width: 630mm

Seat Depth: 460mm

Material: Powder coated round steel tube frame

Features a shaped backrest that gives comfortable lumbar support and a well padded seat that aids easy rising.

