



Mix Medical Cart **Med/Surgical**

The **Med/Surgical Cart** is designed to satisfied with the demanding requirements of virtually every hospital medication distribution system. The Mix Med/Surgical Cart offer the Hospital choices within the innovative Medical Cart line, such as a flexible divider system, removable drawer trays, and a host of organizer accessories nearly arrange the medication and all related supply items. The MIX Med/Surgical Cart also offers the Digital Keyless Security System-I currently.



Digital Keyless Security System-I



Drawer For Tray (only for 3" / 6")

Accessories Options



IV Pole



Tilt Bin Organizers



Sharps Container Set



Glove Dispenser



Waste Container Set(w/ lid)



Bulk Storage Bins

- IV Pole
- Sharps Container Set
- Glove Dispenser
- Waste Container Set(w/ lid)
- Bulk Storage Bins
- Full Locking Narcotics Drawer
- Defibrillator Shelf
- External Power Strip Holder
- Clear Top Mat
- Multi-storage (L)
- Multi-storage (S)
- Bottle Holder
- Cart Handle
- Document Holder

Cart Specifications

Item NO	W		D		H		w/o Casters	
	mm	inch	mm	inch	mm	inch	mm	inch
MX37ANE	690	27-1/8"	518	20-3/8"	1085	45-11/16"	951	37-7/16"
MX34ANE	690	27-1/8"	518	20-3/8"	1006	39-5/6"	872	34-5/16"
MX31ANE	690	27-1/8"	518	20-3/8"	927	36-2/1"	793	31-3/16"
MX28ANE	690	27-1/8"	518	20-3/8"	848	33-8/3"	714	28-1/8"

Drawer Specifications

Item NO	W		D		H	
	mm	inch	mm	inch	mm	inch
A	596	22-3/8"	398	15-11/16"	76	3"
B	596	22-3/8"	398	15-11/16"	155	6"
C	596	22-3/8"	398	15-11/16"	243	9"

ABN: 9400 7000 526

Ph: 03 9878 1399 Fax: 03 9878 7099 National Ph: 1300 65 95 40 National Fax: 1300 65 95 32

Unit 1 Rear 194 – 196 Whitehorse Rd, Blackburn, 3130 VIC